



Founded 1876  
Incorporated 1899

## TOWN OF JEROME, ARIZONA

Jerome Town Hall and Civic Center  
600 Clark Street  
P.O. Box 335, Jerome, AZ 86331  
(928) 634-7943 FAX (928) 634-0715

### APPLICATION FOR DESIGN REVIEW

PLEASE PRINT LEGIBLY SO THAT PROCESSING WON'T BE DELAYED.

APPLICANT'S NAME: JOE WIDMAN

MAILING ADDRESS: P.O. Box 81 JEROME AZ

TELEPHONE: 928-600-0998 EMAIL: JOE@WNF-INC.COM

PROJECT ADDRESS: 146 JUAREZ ST. JEROME AZ

PARCEL NUMBER: \_\_\_\_\_ ZONE DISTRICT: \_\_\_\_\_

APPLICATION FOR (Please describe the project): CUT DOCK ELEVATOR AND  
CONCRETE PAD INSTALLATION FOR HANDICAP ACCESSIBILITY

- I hereby apply for consideration and conditional approval by the Design Review Board.
- I understand that any approval is not valid until application fees are received by the Town.
- I have obtained and reviewed information on the criteria used in evaluation by these bodies and/or reviewed the application provisions in the Town Zoning Ordinance.
- I understand that this application will not be scheduled for consideration until all required materials have been submitted and reviewed.

APPLICANT'S SIGNATURE: Joe Widman DATE: 6-21-18

#### TOWN USE BELOW

RECEIVED FROM: Joe Widman DATE: 6/26/18

Received the sum of \$ 50.00 as: ☐ Check No. \_\_\_\_\_ ☐ Cash ☒ Credit Card

Per Fee schedule - Ordinance 332

BY: Ron Cays

FOR: Charlotte Page

Juarez St.

NE

69'-0"

Water Main

Parking

28'-0"

17'-0"

17'-0"

10'-0"

31'-0"

18'-0"

21'-0"

11'-0"

41'-0"

52'-9 $\frac{5}{8}$ "

SE

NW

120'-0"

Site

48'-0"

50'-0"

SW

Rich St.


N

72'-0"











26








FRONT VIEW



SIDE VIEW



# UNENCLOSED STRAIGHT-THRU PLATFORM (NO PIT)

## VPL-3100/VPL-3200B SERIES VERTICAL PLATFORM LIFT


### TECHNICAL DATA/SPECIFICATIONS

<b>RATED LOAD:</b>	750 lbs maximum.
<b>INPUT POWER SOURCE:</b>	
-AC POWERED UNIT:	110-120 Volt 8 Amp 60 Hz dedicated service
-DC BATTERY POWERED UNIT:	110-120 Volt 3 Amp 60 Hz battery charger.
<b>DRIVE:</b>	
-AC POWERED UNIT:	1 1/2 hp motor, 1750 rpm, 90 VDC, continuous duty.
-DC BATTERY POWERED UNIT:	VPL-3100B: 1/2 hp motor, 1750 rpm, 24 VDC, continuous duty.
	VPL-3200B: 1 hp motor, 1750 rpm, 24 VDC, continuous duty.
	Dual 4L style Poly-V belts and pulleys, 3.94:1 pulley reduction.
<b>INTERMEDIATE REDUCTION:</b>	VPL-3100: 1" dia. ACME screw w/bronze nut and bronze safety back up nut.
<b>FINAL DRIVE:</b>	VPL-3200B: 1.25" dia. ACME screw w/bronze nut and bronze safety back up nut.
<b>MOTOR CONTROLLER:</b>	
-AC POWERED UNIT:	24 VAC Relay control.
-DC BATTERY POWERED UNIT:	24 VDC Relay control.
<b>SPEED:</b>	
-AC POWERED UNIT:	9 feet per minute maximum.
-DC BATTERY POWERED UNIT:	10 feet per minute maximum.

For complete technical specifications please see  
 ILS-00834 "VPL-3100 Series Residential Vertical Platform Lift Technical Specification"  
 ILS-00986 "VPL-3200B Series Residential Vertical Platform Lift Technical Specification"

### PERFORMANCE STANDARDS

USA FOOD & DRUG ADMINISTRATION: CLASS II, 510(K) Exempt, File No. 890.3930,  
 Product Code: PCE  
 ANSI/ASME: A18.1-2014 Section 5 Safety Standards for Platform Lifts and Stairway Chairlifts  
 CSA B355-15 Lifts for Persons with Physical Disabilities (Residential Application)  
 ANSI/ASME: CSA B44.1-14/ASME-A17.5-2014 Elevator and Escalator Electrical Equipment



**ILS-00932**  
**SHEET 1 OF 5**  
**REV. 8 (6440)(8/16/17)(DPG)**

Bruno Independent Living Aids, Inc. 1790 Executive Drive, P.O. Box 84, Oconomowoc, WI 53066  
 Phone (800) 882-8183 Fax (262) 953-5501  
 Email: sales@bruno.com Website: www.bruno.com

**BRUNO**  
 for your independence  
 ©Bruno Independent Living Aids, Inc.

# UNENCLOSED STRAIGHT-THRU PLATFORM (NO PIT)

## VPL-3100/VPL-3200B SERIES VERTICAL PLATFORM LIFT


### ANCHOR POINT LOCATIONS/SLAB DETAIL

#### TECHNICAL SPECIFICATIONS

**FLOOR ATTACHMENT:**  
VPL MUST BE FASTENED TO CONCRETE SLAB USING FOUR (4) 1/2" (3/8" BOLT) X MINIMUM 2-1/2" LONG CONCRETE ANCHORS SUITABLE FOR THE ENVIRONMENT. FOLLOW SELECTED CONCRETE ANCHOR MANUFACTURERS GUIDELINES AND APPLICABLE CODES.

**FLOOR REQUIREMENTS:**  
4" THICK 3500 PSI MINIMUM COMPRESSIVE STRENGTH, REINFORCED CONCRETE SLAB.

For complete technical specifications please see  
ILS-00834 "VPL-3100 Series Residential Vertical Platform Lift Technical Specification"  
ILS-00986 "VPL-3200B Series Residential Vertical Platform Lift Technical Specification"



UNITS WITH STRAIGHT THROUGH PLATFORM WITH PLATFORM GATE


	N	P	Q	R
36" X 48" PLATFORM	7-5/8"	12-1/16"	67-1/2"	54"
36" X 54" PLATFORM	10-5/8"	15-1/16"	73-1/2"	54"
36" X 60" PLATFORM	13-5/8"	18-1/16"	79-1/2"	54"
42" X 60" PLATFORM	13-5/8"	18-1/16"	79-1/2"	60"

**ILS-00932**  
**SHEET 5 OF 5**  
**REV. 8 (6440)(8/16/17)(DPG)**

Bruno Independent Living Aids, Inc. 1780 Executive Drive P O Box 84 Oconomowoc, WI 53066  
Phone (800) 882-8183 Fax (262) 953-5501  
Bruno is a registered trademark of Bruno Independent Living Aids, Inc. All other trademarks are the property of their respective owners.


**BRUNO**  
for your independence  
©Bruno Independent Living Aids, Inc.

# UNENCLOSED STRAIGHT-THRU PLATFORM (NO PIT)



	D	E	F	G	H	J	K	L	M
36"x48" PLATFORM	36-1/4"	33-1/2"	30-3/4"	5-3/4"	49-15/16"	66-1/8"	87-3/16"	31-9/16"	50-13/16"
36"x54" PLATFORM	36-1/4"	33-1/2"	30-3/4"	8-3/4"	55-15/16"	72-1/8"	93-3/16"	31-9/16"	50-13/16"
36"x60" PLATFORM	36-1/4"	33-1/2"	30-3/4"	11-3/4"	63-15/16"	78-1/8"	99-3/16"	31-9/16"	50-13/16"
42"x60" PLATFORM	42-1/4"	39-1/2"	36-3/4"	11-3/4"	61-15/16"	78-1/8"	105-3/16"	34-9/16"	56-13/16"

	A	B	C
VPL-3153	75.9/16"	9.5"	53"
VPL-3175	97.9/16"	11.7"	75"
VPL-32108	148.4/16"	16.5"	123"
VPL-32128	172.1/16"	18.9"	147"
VPL-32148	196.4/16"	21.3"	171"



NOTES: 1) ACCESS RAMP & PLATFORM GATE CAN BE MOUNTED ON EITHER SIDE OF THE PLATFORM  
2) SEE ILS-00938 FOR LANDING GATE DETAIL/ALIGNMENT

**ILS-00932**  
**SHEET 4 OF 5**  
**REV. 8 (6/4/0)(8/16/17)(DPG)**

**Bruno Independent Living Aids, Inc., 1780 Executive Drive, P.O. Box 64, Oconomowoc, WI 53066**  
**Phone (800) 882-8163 Fax (262) 953-5501**

After 18 years, Bruno is now a Division of Johnson & Johnson, Inc. We are pleased to provide the highest quality products and services to our customers. Bruno is now a part of the Johnson & Johnson family.

**BRUNO**  
for your independence  
Bruno Independent Living Aids, Inc.



**President: Art Rios**

**Bus: 928-499-7997**

Email: [arios@cableone.net](mailto:arios@cableone.net)

ROC 202062

**We hearby submit Specifications and Estimates**

We proposee hereby to furnish material and labor-complete in accordance with the above specifications for the sum of:

**Three Thousand six hundred and one**

Dollars

with payments to be made as follows:

### 50% Start 50% Completion

Any alteration from the above specification involving extra costs will be executed

Respectfully

only upon written order, and will become an extra charge over and above the estimate. All

Submitted

## Art Rios

agreements contingent upon strikes, accidents, or delays, beyond our control

Note-this proposal may be withdrawn by us if not accepted within **30** days

## Acceptance of Proposal

The above prices, specifications, and conditions are satisfactory and are hereby accepted. You are authorized to do the work as specified.

**Signature**

Payments will be made as outlined above

**Date of Acceptance**

**Signature**

6/26/18

To whom this may concern,

It is my understanding that Nancy Prochaska would like to install a lift on the north side of her house. This lift would allow her access to her home since she has physical difficulties and can not get to her home via the stairs. As a neighbor I have no problems at all with this installation.

Sincerely,

A handwritten signature in black ink, appearing to read 'Jay Kinsella', with a long horizontal line extending to the right.

Jay Kinsella

40 Rich Street Jerome Az.